



TU Tufts University
CF Core Facility

Peptide Repurification

Contact: _____ Date: _____
 PI Name: _____ Phone: _____
 E-mail: _____ Fax: _____

SHIPPING & BILLING INFORMATION (complete as appropriate):

TUFTS UNIVERSITY	Dept ID:	Project/Grant Number:	
TUFTS MEDICAL CENTER	Cost Center:	Dept:	Box:
ALL OTHER CUSTOMERS			
Shipping:			
Shipping Address:			
Billing Address:			
Purchase Order #:			
Credit Card #: Print Form & Fill-in by hand		EXP:	CVV2:
	<small>*The customer verification number is a non-embossed number located on your card that cannot be obtained easily by using an imprint of your card or simply copying down the account number. Discover/Visa/MasterCard CVV2 numbers are the last three digits printed in the signature field on the back of your card. American Express card CVV2 numbers are four digits long and they are printed above the last four-digit block on the front of the card.</small>		

DO NOT SEND CREDIT CARD INFORMATION BY EMAIL! It is not secure and against credit card company security policy to do so. Please send the completed form by fax or include with samples being sent to us.

FOR YOUR SECURITY, PLEASE SHRED FORM WITH CREDIT INFORMATION AFTER FAXING.



Peptide Repurification

Contact: _____ Date: _____
PI Name: _____ Phone: _____
E-mail: _____ Fax: _____

***** Please include a copy of all the original paperwork that was sent to you.**
(Synthesis Report, Blue Chart Recorder Chromatogram, Mass Spec. etc.)
\$50 Charge if original paperwork is not provided.

1. Peptide Name: _____

Total Number of Purifications: _____ (Maximum of 2 Rounds)

40 mg of Crude Peptide is used for each repurification. Average yield is ~ 5-15 mg of pure peptide.

2. Peptide Name: _____

Total Number of Purifications: _____ (Maximum of 2 Rounds)

40 mg of Crude Peptide is used for each repurification. Average yield is ~ 5-15 mg of pure peptide.

3. Peptide Name: _____

Total Number of Purifications: _____ (Maximum of 2 Rounds)

40 mg of Crude Peptide is used for each repurification. Average yield is ~ 5-15 mg of pure peptide.

4. Peptide Name: _____

Total Number of Purifications: _____ (Maximum of 2 Rounds)

40 mg of Crude Peptide is used for each repurification. Average yield is ~ 5-15 mg of pure peptide.

5. Peptide Name: _____

Total Number of Purifications: _____ (Maximum of 2 Rounds)

40 mg of Crude Peptide is used for each repurification. Average yield is ~ 5-15 mg of pure peptide.

6. Peptide Name: _____

Total Number of Purifications: _____ (Maximum of 2 Rounds)

40 mg of Crude Peptide is used for each repurification. Average yield is ~ 5-15 mg of pure peptide.

Facility use only:

Total Purifications _____ at _____ Shipping _____